

medITECH

MEDICAL GAS EQUIPMENT

Gas Specific Adaptor

Twin Adaptor



- Available with International Fittings.
- 5 Year Warranty.
- Lightweight construction.
- Service Free.

Excellent performance and robustness, coupled with a very lightweight design.

The Meditech Straight Bar Twin Adaptor is designed for use in all Hospitals and Emergency Environments. This device has a Meditech 7 year warranty and is service free for its 10 year recommended life period.

Each Gas Specific Adaptor is available with the probe types and / or the outlet types shown below. If either the probe or outlet type you require is not shown, please contact us to discuss your requirements.

Each Gas Specific Adapter is available with the probe types and / or the outlet types shown below. Other types are available upon request.

PROBE TYPE



BS PROBE



DISS



DIN



AGA



AFNOR



SIS



CZECH



UNIFOR



OUTLET TYPES



BS SCHRADER



DISS



DIN



AGA



AFNOR



SIS



CZECH



UNIFOR




OHIO



SPECIFICATION

Intended Gas:	<ul style="list-style-type: none">• Medical Oxygen (O2) • Medical Air (MA4)• Entonox® (O2/N2O) • Carbon Dioxide (CO2)• Surgical Air (MA7) • Nitrous Oxide (N2O) • Vacuum
Inlet pressure:	4 bar (400 kPa)
Inlet Connections:	Schrader Probe (BS 5682:2015) and other international connections
Outlet Connections (4-BAR):	BS Schrader and Other International Outlets
Intended Life:	10 Years
Service Interval:	Service free
Warranty Period:	7 Years



TO PURCHASE OR FOR FURTHER DETAILS:

Contact us by email or telephone for more details, including specific part numbers, prices, etc.

B.N.O.S. Meditech Ltd.
9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex, CO9 2SZ, England, UK.
TEL: +44 (0) 1787 479 475
FAX: +44 (0) 1787 477747
EMAIL: sales@meditech.uk.com
WEBSITE: www.meditech.uk.com



medITECH

MEDICAL GAS EQUIPMENT

Distributor:

Meditech Gas Specific Adaptors are CE marked in accordance with the European Medical Devices Directive 93/42/EEC and UKCA marked accordance with the UK Medical Devices Regulations 2002. All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice.

B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

Doc 795M Issue C 19/08/2024