



## SUCTION UNIT TROLLEY

The flow-meter suction trolley units allow to perform suction interventions in a rational and effective manner, as they include in one suction system all the necessary devices for a broad range of surgical and medical procedures, providing the operators with the most suitable solution of use in a hospital environment.



This suction unit system allows to connect to the terminal units fitted on the trolley head one or two **EASYVAC**<sup>®</sup> vacuum regulators, preferably completed with the relative **EASYSAFE**<sup>®</sup> safety jar. The slide bracket system on the stand allows to accommodate up to 4 supports for the collection containers of suctioned fluids, disposable or not. The trolley head fits also the catheter container support rings, one for each suction circuit. This to allow the operator to safely store the patient tube after use. The trolley unit fits five antistatic rotating castors, two of which with locking pedal, this to get always an easy, safe and rational positioning. Please refer to the individual product data sheets to see the characteristics of the various devices associated with the trolley.

SIZES (LXWXH)	580x580x870 mm
WEIGHT	5Kg (Vacuum regulators, collection jars and accessories excluded)
TERMINAL UNITS FOR VACUUM REGULATION CONNECTION	2 AFNOR NF-S 90-116 / UNI 9507 / DIN 13260 / BS 5682 / SS 875 24 30 OUTLETS
VACUUM SUPPLY HOSE	EN ISO 5359 YELLOW; LENGTH: 3M
PROBE FOR UNIT CONNECTION TO VACUUM PLANT	AFNOR NF-S 90-116 / UNI 9507 / DIN 13260 / BS 5682 / SS 875 24 30
BRACKET SYSTEM FOR COLLECTION JAR SUPPORT	4 SLIDES WITH DIMENSIONS 25X5 MM

# medITECH

MEDICAL GAS EQUIPMENT

B.N.O.S. Meditech Ltd.  
9 Fifth Avenue, Bluebridge Industrial Estate,  
Halstead, Essex CO9 2SZ, England, UK.  
Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747  
E-mail: sales@meditech.uk.com • www.meditech.uk.com

**TO PURCHASE OR FOR FURTHER DETAILS:**  
**Contact us by email or telephone for more details,**  
**including specific part numbers, prices, etc.**

All items in this range are sold with a return to manufacturer warranty against manufacturing defects. These products are part of a series of medical gas and resuscitation products available from B.N.O.S. Meditech. If the product you require is not shown here, or if you require more information, please contact your distributor or contact us directly. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd Terms and Conditions of sale apply, available on application or from our website.

The Suction Trolley Unit is a registered product, manufactured by Flow Meter Italy.

Distributor: